



# ThunderSTORM: A comprehensive ImageJ plugin for SMLM image processing






Zdeněk Švindrych

Image analysis and data processing in superresolution microscopy workshop 2021  
KONFMI, Prague CZ

August 31, 2021

# Two key concepts: Resolution & Noise

- Resolution





Ernst Abbe (1873):

$$d_{xy} = \frac{\lambda}{2NA}$$

Practically  $\sim 200$  nm


# Two key concepts: Resolution & Noise

- Noise



# Two key concepts: Resolution & Noise

- Resolution and Noise combined



$$\sigma_{xy} = \frac{1}{\sqrt{N}} \sqrt{\sigma_{PSF}^2 + \frac{a^2}{12} + 25 \frac{\sigma_{PSF}^4 b^2}{Na^2}} \approx \frac{\sigma_{PSF}}{\sqrt{N}}$$

“photonic” background plus noise


widefield (Abbe) “resolution”

number of “captured” photons

localization precision (a.k.a. uncertainty)

# Two key concepts: Resolution & Noise

- Resolution and Noise combined




Clarification: Camera offset vs background

$$504 = \text{offset} + \text{readout noise} + \text{background}$$


- Offset = constant value (typically 200 or 500)
  - does not affect the localization precision
- Readout noise = 1 – 2 photoelectrons (sCMOS)
  - < 0.1 photoelectron (EM-CCD)
- Background = any detected photons not originating from the single particle
  - Usually dominant source of noise!

# The Perfect STORM


~ 15000 frames



Particle detection and  
localization




Rendering



Manley *et al*, Current Opinion in Chemical Biology **15**, 2011.

# Not so perfect...



## Critical steps:

- Image filtering
- Detection
- Localization
- Particle filtering

## Bonus:

- Rendering
- Simulation

# SMLM Software

## So, why ThunderSTORM?

- User friendly
- Flexible
- Open
- Published

The SMLM Challenge,  
<http://bigwww.epfl.ch/smlm/challenge2013/index.html?p=participants>


pSMLM-3D	ImageJ	<a href="#">Open access</a>
Auto-Bayes	Stand-alone	<a href="#">Open access</a>
FALCON	Matlab	<a href="#">Open access</a>
MicroManager LM	ImageJ	<a href="#">Open access</a>
SimpleSTORM	Python	<a href="#">Open access</a>
ThunderSTORM	ImageJ	<a href="#">Open access</a>
a-livePALM	Matlab	<a href="#">Access on request</a>
B-recs	C / ImageJ	<a href="#">Open access</a>
Fast-ML-HD	Matlab	
Insight3	Stand-alone	Commercial
L1H	Python	<a href="#">Open access</a>
PeakFit	ImageJ	<a href="#">Open access</a>

SNSMIL	Stand-alone	<a href="#">Open access</a>
SOSplugin	ImageJ	<a href="#">Open access</a>
WaveTracer	Metamorph	Commercial
WTM	Stand-alone	
3D-DAOSTORM	Python / C	<a href="#">Open access</a>
CSSTORM FasterSTORM	Matlab	<a href="#">Open access</a>
GraspJ	ImageJ	<a href="#">Open access</a>
M2LE	ImageJ	<a href="#">Open access</a>
MrSE	Stand-alone	<a href="#">Open access</a>
RapidSTORM	Stand-alone	<a href="#">Open access</a>
DAOSTORM	Python	<a href="#">Open access</a>
PYME	Python	<a href="#">Open access</a>

SFP Estimator (FGPA)	Stand-alone (QT C++)	<a href="#">Open access</a>
GPUgaussMLE gaussMLEv2	Matlab	<a href="#">Open access</a>
MaLiang	ImageJ	<a href="#">Open access</a>
QuickPALM	ImageJ	<a href="#">Open access</a>
SimplePALM	Stand-alone	
Gauss2dcirc	Matlab	<a href="#">Open access</a>
PeakSelector	IDL	
Octane	ImageJ	<a href="#">Open access</a>
Wavelet FluoroBancroft	Matlab	<a href="#">Access on request</a>
Showing 1 to 33 of 33 entries		
© 2018 Biomedical Imaging Group, Ecole Polytechnique Last update: 30 Nov 2018		




# ThunderSTORM



Martin Ovesný



Pavel Křížek





Josef Borkovec



Zdeněk Švindrych



Guy M. Hagen







ThunderSTORM: results

id	fr...	x [nm]	y [nm]	z [nm]	si...	si...	intensi...	offset ...	bkgstd...	chi2	uncert...	detecti...
1	1	972.309	9966...	-6.961	16...	15...	1129.357	17.203	7.085	242962...	9.528	1
2	1	2794...	6583...	-100.38	19...	13...	1517.892	23.068	8.905	383791...	8.968	1
3	1	2897...	7550...	27.032	15...	16...	3578.197	20.912	11.171	604039...	5.083	1
4	1	6564.93	7620...	-353.786	33...	15...	2199.385	25.737	8.495	349311...	11.146	1
5	1	7082...	2957...	-63.09	18...	14...	1414.138	18.272	7.866	299472...	8.618	1
6	1	7081...	6146...	-226.521	25...	13...	2896.539	38.236	11.073	593394...	7.472	1
7	1	7503.76	13196...	-71.456	18...	14...	1717.186	20.168	12.612	769921...	9.825	1
8	1	8653.55	12139...	-157.855	21...	13...	1984.246	16.686	8.091	316883...	7.481	1
9	1	9202...	1807...	187.182	14...	22...	1957.3	20.858	8.338	336475...	8.176	1
10	1	10091...	3430...	143.474	14...	20...	1350.3	22.218	7.359	262128...	9.461	1
11	1	10263...	7047...	27.624	15...	16...	3445.96	30.986	11.884	683549...	5.377	1
12	1	10546...	645.07	0.702	16...	15...	2275.419	30.799	15.253	112598...	8.625	1
13	1	2866...	4961...	-78.758	18...	13...	6845.859	25.493	20.687	4.365	2	
14	1	10001...	4446...	54.982	15...	16...	6341.438	22.522	27.104	5.552	3	
15	2	2156...	2744...	-118.054	19...	13...	1117.588	29.296	7.726	288873...	10.715	1
16	2	2270...	10662...	-159.182	21...	13...	2021.478	19.81	8.577	356025...	7.601	1
17	2	7553...	13858...	-104.47	19...	13...	988.367	13.635	6.287	191290...	10.405	1
18	2	9520...	6632...	-104.694	19...	13...	1232.662	26.11	10.669	550922...	11.94	1

Filter: Density filter Remove duplicates Merging Drift correction Z-stage offset

Post-processing history: -

Buttons: Preview, Defaults, Plot histogram, Visualization, Import, Export



# ThunderSTORM paper

## ThunderSTORM: a comprehensive ImageJ plug-in for PALM and STORM data analysis and super-resolution imaging

Martin Ovesný, Pavel Krížek, Josef Borkovec, Zdeněk Švindrych and Guy M. Hagen\*

Institute of Cellular Biology and Pathology, First Faculty of Medicine, Charles University in Prague, Prague 12800, Czech Republic

Associate Editor: Jonathan Wain

### ABSTRACT

**Summary:** ThunderSTORM is an open-source, interactive and modular plug-in for ImageJ designed for automated processing, analysis and visualization of data acquired by single-molecule localization microscopy methods such as photo-activated localization microscopy and stochastic optical reconstruction microscopy. ThunderSTORM offers an extensive collection of processing and post-processing methods so that users can easily adapt the process of analysis to their data. ThunderSTORM also offers a set of tools for creation of simulated data and quantitative performance evaluation of localization algorithms using Monte Carlo simulations.

**Availability and implementation:** ThunderSTORM and the online documentation are both freely accessible at <https://code.google.com/p/thunder-storm/>.

**Contact:** guy.hagen@ffl.cuni.cz

**Supplementary information:** Supplementary data are available at Bioinformatics online.

Received on February 4, 2014; revised on April 9, 2014; accepted on April 13, 2014

### 1 INTRODUCTION

Single-molecule localization microscopy (SMLM) methods such as stochastic optical reconstruction microscopy (STORM; Rust *et al.*, 2006) and photo-activated localization microscopy (PALM; Betzig *et al.*, 2006) have recently emerged to overcome the diffraction barrier, offering ~10 times higher lateral resolution and the possibility of 3D imaging by various approaches. In SMLM, a super-resolution image is reconstructed from a sequence of diffraction-limited images of sparsely distributed single photoswitchable molecules. As the sequence is usually long (thousands of images) and the positions of the molecules have to be estimated systematically with sub-diffraction precision, it is crucial to use specialized software for processing the data.

We present ThunderSTORM, an open-source, interactive, modular and platform-independent software, which provides a complete set of tools for automated processing, analysis and visualization of data acquired by SMLM methods. Our philosophy in developing ThunderSTORM has been to offer an extensive collection of processing and post-processing methods, which were developed based on extensive testing with both real and

simulated data. We also provide a detailed description of the implemented methods and algorithms (Supplementary Note), as well as a user's guide.

### 2 FEATURES AND METHODS

Most software tools currently available for SMLM data processing typically use only one particular algorithm for detection and localization of molecules. ThunderSTORM offers many different processing and post-processing methods so that users can adapt the analysis to their data. This approach can lead to higher quality results than existing solutions. Experienced users may use any combination of the available methods; however, we have designed the software's default settings to produce good results on many of the datasets we have experimented with.


#### 2.1 Raw data processing

Approximate molecular positions can be determined, in combination with a variety of feature-enhancing low-pass and band-pass filters (Krížek *et al.*, 2011; Izeddin *et al.*, 2012), by detection of local maxima, non-maximum suppression or calculation of the centroid of connected components of segmented objects. A feature exclusively unique to ThunderSTORM is the possibility of specifying the threshold for detection of molecules using a mathematical expression with quantities based on raw or filtered images. This allows computing the threshold value systematically for unknown input images with, for example, low signal to noise ratio, or where the global intensity slowly fluctuates. ThunderSTORM also offers a preview function to help visualize the detected molecules with the chosen combination of data processing settings.

Sub-diffraction localization of molecules is accomplished by computing the centroid of local neighborhoods, by a radial symmetry approach (Parthasarathy, 2012), or by fitting a suitable PSF model using standard or weighted non-linear least-squares methods, or using maximum-likelihood estimation (Mortensen *et al.*, 2010). Users may also choose not to use any of the methods, thereby using the approximate localizations from the previous step. The uncertainty of the localization of molecules is calculated according to Thompson *et al.* (2002), or according to Quan *et al.* (2010) if EMCCD cameras are used.

Super-resolution 3D imaging is accomplished by an astigmatism approach (Huang *et al.*, 2008). An integral part of this feature is the software's calibration tool, in which a Z-stack of astigmatic images of sub-diffraction fluorescent beads is used to establish parameters for determining the axial position of each molecule.

Efforts to accelerate the acquisition process in SMLM have involved increasing the density of photoactivated fluorophores. In this case, ThunderSTORM uses an algorithm based on fitting of multiple emitters (Huang *et al.*, 2011).



**Fig. 1.** Simulations and SMLM reconstruction with ThunderSTORM. (a) Example of a mask used for generating simulated SMLM data. The gray-scale intensity values are interpreted as molecular densities within a user-specified range. (b) SMLM reconstruction of a simulated dataset. (c) Detail of (b). (d) Widefield fluorescence image of an A431 epidermoid carcinoma cell expressing the membrane protein mChitrine-erbB3. (e) SMLM reconstruction using the default settings. (f) Detail of (e). SMLM imaging was performed as previously described (Krížek *et al.*, 2011).

#### 2.2 Post-processing and visualization

Post-processing routines offered by ThunderSTORM can eliminate molecules with poor localization or other user-defined criteria, merge molecules reappearing in subsequent frames, remove duplicated molecules obtained in multiple emitter analysis (Huang *et al.*, 2011), correct molecular positions for lateral drift of the sample using fiducial markers or using cross-correlation methods (Modziński *et al.*, 2011) and correct the absolute axial position of the molecules when the data were acquired in multiple Z-stack positions (Huang *et al.*, 2008). Users can also select a region of interest to export only the localized molecules and their parameters from the region. Post-processing includes a live preview.

Visualization involves creation of a new high-resolution image based on the previously obtained sub-diffraction molecular coordinates. Several methods have been implemented for visualization such as Gaussian rendering and a 2D histogram with an option of jittering (Krížek *et al.*, 2011). ThunderSTORM also introduces a new visualization method based on an average shifted histogram approach (Scott, 1985). This method provides similar results as the Gaussian rendering, but is orders of magnitude faster.

#### 2.3 Simulation engine and performance evaluation

ThunderSTORM is capable of generating realistic sequences of SMLM-like images in which the ground-truth positions of the molecules are known. A grayscale mask can be used to vary the spatial density of molecules [Fig. 1(a–c)]. When the localization data and the ground-truth positions of molecules are available, ThunderSTORM can quantitatively evaluate the performance of localization algorithms (see Supplementary Note Sections 8 and 9). This allows users to perform sophisticated Monte Carlo simulations (Krížek *et al.*, 2011) (see User's Guide Sections 8–10).

### 3 SUMMARY

ThunderSTORM introduces several new features and concepts for 2D and 3D SMLM data analysis. The software combines several algorithms for SMLM analysis into one comprehensive environment. One of the main features is the ability to process the data using any combination of the implemented feature-enhancing, spot detection and fitting methods. An important feature in ThunderSTORM is the possibility of specifying the threshold for detection of molecules using mathematical expressions. This allows users to systematically maximize the efficiency of molecule detection in the raw data by searching for the optimum combination, which may vary from experiment to experiment. ThunderSTORM also offers a much higher degree of user interactivity during data post-processing compared with other SMLM software packages, and introduces a new and fast visualization method that creates high-quality results. A realistic data generator within ThunderSTORM allows users to run multidimensional Monte Carlo simulations to evaluate the performance of localization methods. We have found ThunderSTORM's flexibility and performance to be of critical importance when analyzing data with low molecular brightness, which we encountered when imaging A431 cells expressing mChitrine-erbB3 (Krížek *et al.*, 2011) [Fig. 1(d–f)].

**Funding:** This work was supported by the Czech Science Foundation [P304/09/1047, P205/12/P392, P302/12/G157, 14-15272P; Charles University [Prvok/ILF/1, UNCE 204022]; the European Regional Development Fund [OPVK CZ.2.16/3.1.00/24/010, BIOCEV CZ.1.05/1.1.00/02.0109]; and the European Social Fund [OPVK CZ.1.07/2.3.00/30.0030].

**Conflict of Interest:** none declared.

### REFERENCES


- Betzig, E. *et al.* (2006) Imaging intracellular fluorescent proteins at nanometer resolution. *Science*, **313**, 1480–1485.
- Huang, B. *et al.* (2008) Three-dimensional super-resolution imaging by stochastic optical reconstruction microscopy. *Science*, **319**, 810–813.
- Huang, B. *et al.* (2011) Simultaneous multiple-emitter fitting for single molecule super-resolution imaging. *Bioimed. Opt. Express*, **2**, 1377–1393.
- Izeddin, I. *et al.* (2012) Wavelet analysis for single molecule localization microscopy. *Opt. Express*, **20**, 2881–2895.
- Krížek, P. *et al.* (2011) Minimizing detection errors in single molecule localization microscopy. *Opt. Express*, **19**, 3226–3235.
- Modziński, M.J. *et al.* (2011) Sample drift correction in 3D fluorescence photo-activation localization microscopy. *Opt. Express*, **19**, 15009–15019.
- Mortensen, K.J. *et al.* (2010) Optimal localization analysis for single-molecule tracking and super-resolution microscopy. *Nat. Methods*, **7**, 377–381.
- Parthasarathy, R. (2012) Rapid, accurate particle tracking by calculation of radial symmetry centers. *Nat. Methods*, **9**, 724–725.
- Quan, T. *et al.* (2010) Localization capability and limitation of electron-multiplying charge-coupled, scientific complementary metal-oxide semiconductor, and charge-coupled devices for super-resolution imaging. *J. Biomed. Opt.*, **15**, 066005.
- Rust, M.J. *et al.* (2006) Sub-diffraction-limit imaging by stochastic optical reconstruction microscopy (STORM). *Nat. Methods*, **3**, 793–795.
- Scott, D.W. (1985) Averaged shifted histograms: effective nonparametric density estimators in several dimensions. *Ann. Stat.*, **13**, 1034–1040.
- Thompson, R.E. *et al.* (2002) Precise nanometer localization analysis for individual fluorescent probes. *Biophys. J.*, **82**, 2775–2783.

Comes with:

- User Guide – 27 pages
- Methods and algorithms – 40 pages
- Online resources <https://github.com/zitmen/thunderstorm/wiki>
- Hundreds of thousands of lines of Java code !! (most of them are just “{” and “}” ...)
- Martin's PhD Thesis – 150 pages!

\*To whom correspondence should be addressed.


# Data processing pipeline



# Image acquisition



Won't go into any detail here.

- STORM microscope (strong lasers, sensitive camera)
- Blinking sample!!!
- Minimize background
- Proper sampling (Nyquist) – 100 nm pixels or smaller
- High labeling density (achievable for 1D structures or small clusters)
- 3D – astigmatic lens,  
z-range can be extended with a z-stack



# Open image stack

- Handled by ImageJ (FIJI, BioFormats)
- Note, native .tif is 5x faster than BioFormats OME TIFF
- Virtual stack supported!
- Check pixel size
- Multichannel images are processed separately





# Camera setup

- Basic camera info for correct noise handling
- Small effect on localization results (MLE)
- Strong effect on the reported localization precision
- Strong effect on the reported intensities of individual molecules

Pixel size [nm]:	117.0
Photoelectrons per A/D count:	5.0
Base level [A/D counts]:	400.0
<input checked="" type="checkbox"/> EM gain:	100.0

Defaults ? OK Cancel

Value of totally dark pixels  
a.k.a. Offset

# Run analysis

- The main dialog...

Plugins → ThunderSTORM → Run analysis

The screenshot shows the 'Run analysis' dialog box in ThunderSTORM. It is divided into several sections: 'Camera' (with a 'Camera setup' button), 'Image filtering' (with a 'Filter' dropdown set to 'Wavelet filter (B-Spline)', 'B-Spline order' set to 3, and 'B-Spline scale' set to 2.0), 'Approximate localization of molecules' (with a 'Method' dropdown set to 'Local maximum', 'Peak intensity threshold' set to 'std(Wave.F1)', and 'Connectivity' set to '8-neighbourhood'), 'Sub-pixel localization of molecules' (with a 'Method' dropdown set to 'PSF: Integrated Gaussian', 'Fitting radius [px]' set to 3, 'Fitting method' set to 'Least squares', and 'Initial sigma [px]' set to 1.6), 'Multi-emitter fitting analysis' (with a 'Multi-emitter fitting analysis' checkbox set to 'enable', 'Maximum of molecules per fitting region' set to 5, 'Model selection threshold (p-value)' set to 1.0E-6, and checkboxes for 'Same intensity for all molecules' and 'Limit intensity range [photons]' set to 500:2500), and 'Visualisation of the results' (with a 'Method' dropdown set to 'Averaged shifted histograms', 'Magnification' set to 5.0, 'Update frequency [frames]' set to 50, '3D' checkbox set to 'off', 'Colorize z-stack' checkbox set to 'off', 'Z range (from:step:to) [nm]' set to -500:100:500, 'Lateral shifts' set to 2, and 'Axial shifts' set to 2). At the bottom are buttons for 'Defaults', 'Preview', 'Ok', and 'Cancel'. Red arrows point to various sections: 'Camera parameters setup.' points to the 'Camera' section; 'Image filtering and feature enhancement.' points to the 'Image filtering' section; 'Finding approximate position of molecules.' points to the 'Approximate localization of molecules' section; 'Setting the threshold.' points to the 'Peak intensity threshold' field; 'Sub-pixel 2D/3D localization of molecules.' points to the 'Sub-pixel localization of molecules' section; 'Crowded-field problem.' points to the 'Multi-emitter fitting analysis' section; 'Instant visualization of the super-resolution image during data analysis.' points to the 'Visualisation of the results' section; and 'Preview of molecular localizations with current settings.' points to the 'Preview' button. A green arrow points to the 'Multi-emitter fitting analysis' section.

Camera parameters setup.

Image filtering and feature enhancement.

Finding approximate position of molecules.

Setting the threshold.

Sub-pixel 2D/3D localization of molecules.

Crowded-field problem.

Instant visualization of the super-resolution image during data analysis.

Preview of molecular localizations with current settings.


# 1. Image filtering

Remove noise and background without affecting the image resolution (uh?)

- Averaging filter (square box)
- Gaussian filter
- Lowered Gaussian (used in DAOSTORM)
- Difference of Gaussians
- Median
- Wavelet (recommended)

Note:

The filter is only used to find candidate molecules.  
The subsequent fitting routines operate on the raw data!






## 2. Local maxima search

Approximate positions of the candidate molecules

- Three methods available based on published papers
- Results are quite similar. Local maximum with 8-connected neighborhood shows best results in tests...
- Threshold – single value for the whole frame.
  - actual threshold value reported in Results (ImageJ)
  - flexible formula parser
  - `std(Wave.F1)*2` (noisy data)
  - `std(Wave.F1)*4` (clean data)



### Built-in operators

- `a + b`
- `a - b`
- `a * b`
- `a / b`
- `a % b`
- `a ^ b`

### Built-in functions

- `var(x)`
- `std(x)`
- `mean(x)`
- `median(x)`
- `min(x)`
- `max(x)`
- `sum(x)`
- `abs(x)`

### Variables provided by different feature enhancement filters

All the filters provide these variables:

- `I` – input image without any changes
- `F` – final filtered image


Variables representing the filters:

- `Box.F` – Averaging (Box) filter
- `Med.F` – Median filter
- `Gauss.F` – Gaussian filter
- `LowGauss.F` – Lowered Gaussian filter
- `DoG.F` – Difference-of-Gaussians filter
  - `DoG.G1` – Input image filtered by 1<sup>st</sup> Gaussian function of DoG filter
  - `DoG.G2` – Input image filtered by 2<sup>nd</sup> Gaussian function of DoG filter
- `DoB.F` – Difference of averaging filters
  - `DoB.B1` – Input image filtered by 1<sup>st</sup> Box of DoB filter
  - `DoB.B2` – Input image filtered by 2<sup>nd</sup> Box of DoB filter
- `Wave.F` – Wavelet filter
  - `Wave.F1` – 1<sup>st</sup> Wavelet level of the input image
  - `Wave.F2` – 2<sup>nd</sup> Wavelet level of the input image

# 3. Subpixel localization

## The core of ThunderSTORM

- Fitting methods: Gaussian, Integrated Gaussian, Elliptical Gaussian (3D, needs 3D calibration)
  - combined with Least Squares, Weighted-LS, or Maximum Likelihood
- Non-fitting methods: Centroid, Radial symmetry
- Multi-emitter fitting
  - computationally intensive ( = slow!)
  - parameters to limit the number of particles/region
- Localized particles are rendered in a superresolution image on the fly (more on Visualization later)



# The Results

Sortable columns

ThunderSTORM: results

id	frame	x [nm]	y [nm]	z [nm]	sigma1 [nm]	sigma2 [nm]	intensity [photon]	offset [photon]	bkgstd [photon]	chi2	uncertainty [nm]	detections
164253	13545	11668.133	13985.968	-116.226	199.19	136.337	2708.894	17.909	9.925	4767579.916	6.217	1
164254	13545	14013.344	12750.633	164.394	143.347	210.813	1954.007	16.029	7.349	2614093.167	7.432	1
164255	13545	14279.572	2473.196	138.473	143.951	199.348	1390.258	20.769	12.08	7062292.89	12.5	1
164256	13545	5541.074	10387.029	-154.819	215.466	134.75	4925.082	24.829	22.838	◆	6.731	2
164257	13545	7899.836	5743.203	88.588	147.738	179.925	7391.769	29.154	27.865	◆	5.078	3
164258	13545	9527.411	11012.639	-99.692	192.904	137.708	4845.302	20.425	20.304	◆	5.827	2
164259	13546	978.339	9202.801	-273.148	277.778	142.384	1943.128	28.482	13.686	9066230.981	13.557	1
164260	13546	3088.897	3948.948	-271.685	276.889	142.17	1991.757	48.529	22.213	23881409.49	19.849	1
164261	13546	4125.555	11177.337	33.542	155.869	162.472	1302.911	19.755	7.613	2805272.351	8.941	1
164262	13546	6474.712	9928.363	-51.185	176.575	143.859	1686.989	23.869	10.953	5806254.382	8.985	1
164263	13546	6803.453	10993.918	-79.363	185.656	139.878	1128.19	19.539	6.776	2222228.717	9.577	1
164264	13546	7581.077	1541.244	349.093	165.875	319.532	4237.493	25.854	19.89	19147386.572	11.786	1
164265	13546	8919.125	8509.081	-201.911	237.892	135.399	2072.828	31.027	8.519	3512809.184	7.988	1
164266	13546	11658.317	13983.667	-136.465	207.418	135.195	2201.414	13.445	8.475	3476107.93	6.884	1
164267	13546	13527.045	5654.882	157.187	143.422	207.531	1007.086	21.842	6.793	2233564.943	11.587	1
164268	13546	1541.363	11467.352	-87.518	188.499	138.942	5172.26	19.897	18.081	◆	5.073	2
164269	13546	7581.838	1527.354	353.678	167.057	322.858	13112.9	25.93	60.864	◆	11.108	3
164270	13546	9767.863	10032.998	214.787	144.877	235.826	9513.656	30.364	35.215	◆	6.055	2

Filter Density filter Remove duplicates Merging Drift correction Z-stage offset

Filter:  Apply Restrict to ROI

Post-processing history: - Reset

☒ Preview Defaults Plot histogram Visualization Import Export


Filtering and post-processing

Undo history

Help!!!

# Filtering results


## Handy Histograms



# Filtering results

## Density filter

- Removes spurious localizations far from other localizations (when the density of localizations is very low)
- Sort of noise removal




The screenshot shows a software interface with a tabbed menu at the top containing 'Filter', 'Density filter', 'Remove duplicates', 'Merging', 'Drift correction', and 'Z-stage offset'. The 'Density filter' tab is active. Below the tabs, there are two input fields: 'Distance radius [nm]:' with the value '50.0' and 'Minimum number of neighbors in the radius:' with the value '5'. To the right of these fields are two radio buttons, '2D' (which is selected) and '3D'. Further right is a blue question mark icon. Below the input fields are two buttons: 'Apply' and 'Reset'. At the bottom left, there is a label 'Post-processing history: -'.

# Filtering results

## Remove duplicates

- Molecules with the same coordinates in a given frame are removed
- Duplicates may arise due to overlap of the fitting regions
- Usually not a dramatic effect...



The screenshot shows a software interface with a tabbed menu at the top containing 'Filter', 'Density filter', 'Remove duplicates', 'Merging', 'Drift correction', and 'Z-stage offset'. The 'Remove duplicates' tab is active. Below the tabs, the text reads 'Remove molecules that converged to the same position. Distance threshold [nm]:' followed by a text input field containing 'uncertainty' and a blue question mark icon. To the right of the input field is an 'Apply' button. Below this, the 'Post-processing history:' section shows 'Remove duplicates' with a checked checkbox. To the right of this section is a 'Reset' button. At the bottom left, there is a 'Preview' checkbox which is checked. To the right of the 'Preview' checkbox are five buttons: 'Defaults', 'Plot histogram', 'Visualization', 'Import', and 'Export'. At the very bottom of the dialog, a status bar displays the text '557 out of 276711 items were filtered out'.

# Filtering results

## Merging of molecules

- Molecules at the same position in subsequent frames are merged
- Molecule “on” time is longer than camera exposure time

The screenshot shows a software window with several tabs: Filter, Density filter, Remove duplicates, Merging (selected), Drift correction, and Z-stage offset. In the Merging tab, there are two input fields: 'Maximum distance [units of x,y]:' with the value 20, and 'Max. frames per molecule (0 = unlimited):' with the value 0. Below these is another input field 'Maximum off frames:' with the value 1. To the right of these fields is a 'Merge' button. Below the input fields is a 'Post-processing history:' section showing 'Merging' with a checked checkbox. To the right of this is a 'Reset' button. At the bottom of the window, there is a status bar that reads '276711 molecules were merged into 214702 molecules'. Above the status bar are several buttons: 'Preview' (checked), 'Defaults', 'Plot histogram', 'Visualization', 'Import', and 'Export'.

Filter Density filter Remove duplicates Merging Drift correction Z-stage offset

Maximum distance [units of x,y]: 20 Max. frames per molecule (0 = unlimited): 0

Maximum off frames: 1

Post-processing history: Merging ☒




☒ Preview Defaults Plot histogram Visualization Import Export

276711 molecules were merged into 214702 molecules

# Filtering results

## Drift correction

- Cross-correlation – usually quite reliable
- Fiduciary markers – less tested






# Filtering results

## Z-stage offset

- Extends the useful imaging depth

Filter	Density filter	Remove duplicates	Merging	Drift correction	Z-stage offset
Frames per one Z-stage position: 100		Z-stage step [nm]: 500		?	
Number of Z-stage positions: 4		Initial Z-stage offset [nm]: 0		Apply	
Post-processing history:					Reset
<input checked="" type="checkbox"/> Preview		Defaults	Plot histogram	Visualization	Import Export
Drift correction applied.					




# Visualization

Renders the superresolution image from the localization coordinates in the ThunderSTORM results table

Several methods available:


- Normalized Gaussian – draws a Gaussian in place of the molecule
- Shifted histograms – quicker approximation to Gauss
- Histogram – 100 jittered 1x1 pixel dots to the image
- Scatterplot – adds 1x1 pixel dots to the image

3D image rendered as stack, can be colored



# Visualization



3D projection in ImageJ



# 3D STORM

Z-localization requires modified PSF

- ThunderSTORM currently supports astigmatic PSF (cylindrical lens)
- Calibration performed on beads
  - select higher thresholds



Parameters of the fitted calibration curves and the angle of the astigmatic lens with respect to the camera.


Calibration z-stack, where beads used for calibration are marked with red circles and slice-by-slice sub-pixel localizations are marked with blue crosses.

Fitted calibration curves.

# Other features


## Simulation engine

- Generates SMLM time series images based on grayscale masks and desired blinking statistics
- Useful for evaluation of processing parameters (we have the ground truth)



## Coordinate-based colocalization



- Evaluating relative coordinates of molecules from two-channel images
- Also useful in cluster analysis




AF647-EDU vs FU (anti-BrDU)  
Sample courtesy Eugene Smirnov



# Some resulting images



# Comparison with other software




The SMLM Challenge,  
<http://bigwww.epfl.ch/smlm/challenge2013/>



# Conclusion

- ThunderSTORM is very capable, easy to use software for SMLM data processing
- Open source, verified algorithms
- Freely adjustable parameters to suit most SMLM data
- No GPU acceleration
- Currently only Astigmatism-based 3D SMLM
- Not being actively developed (at least not by us)...



Dual Objective STORM, U2OS cell  
courtesy Martin Ovesný



# Credits

Guy M. Hagen

Martin Ovesný

Pavel Křížek

Josef Borkovec

Eugene Smirnov

Ivan Raška

The SMLM Challenge organizers at <http://bigwww.epfl.ch/smlm/challenge2013/>